



Infinity Buccal Tubes

First Molar non-convertible



Convertible



Features & Benefits


- Bondable, single, non-convertible or convertible tubes which are cast with tooth ID.
- Two-Piece Tubes with a mesh base ensures exceptional bonding strength.
- Low-Profile provides ultimate patient comfort and second molar tubes are smaller in size for partially erupted cases.
- Indent in the base makes mesio-distal placement easier.
- Available in .018" Roth & .022" M.B.T / all other prescriptions, double and triple tubes available on request.

System	Tooth	Torque	Offset	In/out	Width
Roth	UR	-14°	10°	0.4	4.2
	UL	-14°	10°	0.4	4.2
	LL	-25°	10°	0.4	4.2
	LR	-25°	10°	0.4	4.2
MBT	UR	-14°	10°	0.4	4.2
	UL	-14°	10°	0.4	4.2
	LL	-20°	10°	0.4	4.2
	LR	-20°	10°	0.4	4.2
Edgewise	UR/LL	0°	10°	0.4	4.2
	UL/LR	0°	10°	0.4	4.2



Infinity Buccal Tubes

Second Molar non-convertible



Features & Benefits

- Bondable, single, non-convertible or convertible tubes which are cast with tooth ID.
- Two-Piece Tubes with a mesh base ensures exceptional bonding strength.
- Low-Profile provides ultimate patient comfort and second molar tubes are smaller in size for partially erupted cases.
- Indent in the base makes mesio-distal placement easier.
- Available in .018" Roth & .022" M.B.T / all other prescriptions, double and triple tubes available on request.
- Minimum quantity order of 5 tubes per quadrant.

System	Tooth	Torque	Offset	In/out	Width
Roth	UR	-14°	10°	0.4	3.2
	UL	-14°	10°	0.4	3.2
	LL	-25°	4°	0.4	3.2
	LR	-25°	4°	0.4	3.2
MBT	UR	-14°	10°	0.4	3.2
	UL	-14°	10°	0.4	3.2
	LL	-10°	0°	0.4	3.2
	LR	-10°	0°	0.4	3.2