

Smile Safari™ Intraoral Elastics Chart

	Light - 2.5 oz.	Medium - 4.5 oz.	Heavy - 6.5 oz.
1/8"	 Giraffe 1/8" (3.2mm) Light 2.5 oz.	 Chimpanzee 1/8" (3.2mm) Medium 4.5 oz.	 Python 1/8" (3.2mm) Heavy 6.5 oz.
3/16"	 Toucan 3/16" (4.8mm) Light 2.5 oz.	 Zebra 3/16" (4.8mm) Medium 4.5 oz.	 Python 3/16" (4.8mm) Heavy 6.5 oz.
1/4"	 Gazelle 1/4" (6.4mm) Light 2.5 oz.	 Lion 1/4" (6.4mm) Medium 4.5 oz.	 Rhino 1/4" (6.4mm) Heavy 6.5 oz.
5/16"	 Leopard 5/16" (7.9mm) Light 2.5 oz.	 Gorilla 5/16" (7.9mm) Medium 4.5 oz.	 Rhino 5/16" (7.9mm) Heavy 6.5 oz.
3/8"	 Vulture 3/8" (9.4mm) Light 2.5 oz.	 Warthog 3/8" (9.4mm) Medium 4.5 oz.	 Crocodile 3/8" (9.4mm) Heavy 6.5 oz.