

BUCCAL TUBES

First Molar Bonding Buccal Tubes

Single Tubes



System	Tooth	Torque	Offset	In/out	Width
Roth	UR	-14°	10°	0.4	4.2
	UL	-14°	10°	0.4	4.2
	LL	-25°	10°	0.4	4.2
	LR	-25°	10°	0.4	4.2
MBT	UR	-14°	10°	0.4	4.2
	UL	-14°	10°	0.4	4.2
	LL	-20°	10°	0.4	4.2
	LR	-20°	10°	0.4	4.2
Edgewise	UR/LL	0°	10°	0.4	4.2
	UL/LR	0°	10°	0.4	4.2

Available in .018" or .022" Roth or MBT, convertible and double or triple tubes.



Second Molar Bonding Buccal Tubes

Single Tube

System	Tooth	Torque	Offset	In/out	Width
Roth	UR	-14°	10°	0.4	3.2
	UL	-14°	10°	0.4	3.2
	LL	-25°	4°	0.4	3.2
	LR	-25°	4°	0.4	3.2
MBT	UR	-14°	10°	0.4	3.2
	UL	-14°	10°	0.4	3.2
	LL	-10°	0°	0.4	3.2
	LR	-10°	0°	0.4	3.2

Available in .018" or .022" Roth or MBT.

Please note: Minimum quantity order of 5 tubes per quadrant