




Clear Self-Ligating Brackets

 **Diamond Clear**TM
SELF-LIGATING BRACKETS





Clear Self-Ligating Brackets

Features & Benefits

- The Diamond Clear Self-Ligating (SL) Bracket is the latest addition to our aesthetic range and is manufactured using the latest unique laser technique.
- The Diamond-like Clarity ensures optimal aesthetics.
- Unique design of sliding gate prevents accidental pull-out.
- The Diamond Clear SL Bracket is manufactured from a very innovative and new Premium Material which is the next generation material beyond polycrystalline alumina. The Diamond Clear SL Bracket is unbelievably clear and does not stain like the common polyurethane and polycarbonate materials and does not abrade opposing teeth like ceramic brackets.
- The material of the Diamond Clear SL Brackets is hard enough not to deform during treatment like standard composite brackets tend to do.
- Due to the strength of the material, bracket fractures during debonding are minimised.
- The bracket also features a mechanical locking fishgrate base. There is no additional coating or silane needed.
- The special base provides excellent shear forces and easy debonding.

Special Base for easy bonding and debonding



FRICTION FORCES

Due to the new material composition the brackets show outstanding ultra-low frictional forces for premium sliding mechanics compared to all other used bracket materials on the market. This allows full control of the treatment forces.

