

Self-Ligating Brackets

Diamond QMTM

SELF-LIGATING BRACKETS




Features & Benefits

- Torque-in-Base design
- Uses passive ligating system
- Quick and easy wire changes
- Optimised slot depth, which is closer to the tooth surface
- Use the Diamond QM optimally with CU Alloy NiTi

Diamond QM Torque

	Tooth	Torgue°	Tip°	Offset°	In-Out/mm	Wide/mm	High/mm
MAXILLARY	U1	17	4	0	0.71	3.0	2.3
	U2	10	8	0	0.97	2.8	2.5
	U3	-7	8	0	0.36	3.0	2.0
	U3 Hook	-7	8	0	0.36	3.0	2.0
	U4*5	-7	0	0	0.36	3.0	2.0
	U4*5 Hook	-7	0	0	0.36	3.0	2.0
MANDIBULAR	L1*2	-6	0	0	1.07	2.6	2.7
	L3	-6	3	0	0.46	3.0	2.0
	L3 Hook	-6	3	0	0.46	3.0	2.0
	L4	-12	2	0	0.28	3.0	1.8
	L4 Hook	-12	2	0	0.28	3.0	1.8
	L5	-17	2	0	0.28	3.0	1.8
	L5 Hook	-17	2	0	0.28	3.0	1.8



ORTHOSHOP

Your first step to the perfect smile