


# Super Mini Twin Brackets



## Features & Benefits

- Comes with a 80 Gauge mesh base (rather than a cast MIM base) which has also been sandblasted thereby further increasing retention.
- They have a chamfered slot and perfectly round edges for patient comfort.
- Manufactured from the best quality 17-4 Stainless Steel and MIM (Metal Injection Mould) technology.
- Available in mini size.
- Available in .018" / .022" Roth & M.B.T (see Torque and Angulation information below), other prescriptions and sizes available on request

### MBT Bracket Slot

| Tooth              | Torque | Angulation |
|--------------------|--------|------------|
| UL1 / 0.022/0.018  | +17°   | 4°         |
| UR1 / 0.022/0.018  | +17°   | 4°         |
| UL2 / 0.022/0.018  | +10°   | 8°         |
| UR2 / 0.022/0.018  | +10°   | 8°         |
| UL3 / 0.022/0.018  | 0°     | 8°         |
| UR3 / 0.022/0.018  | 0°     | 8°         |
| UL45 / 0.022/0.018 | -7°    | 0°         |
| UR45 / 0.022/0.018 | -7°    | 0°         |
| LL12 / 0.022/0.018 | -6°    | 0°         |
| LR12 / 0.022/0.018 | -6°    | 0°         |
| LL3 / 0.022/0.018  | -6°    | 3°         |
| LR3 / 0.022/0.018  | -6°    | 3°         |
| LL4 / 0.022/0.018  | -12°   | 0°         |
| LR4 / 0.022/0.018  | -12°   | 0°         |
| LL5 / 0.022/0.018  | -17°   | 0°         |
| LR5 / 0.022/0.018  | -17°   | 0°         |

TOLERANCE ±0.5

### Roth Bracket Slot

| Tooth              | Torque | Angulation |
|--------------------|--------|------------|
| UL1 / 0.022/0.018  | +12°   | 5°         |
| UR1 / 0.022/0.018  | +12°   | 5°         |
| UL2 / 0.022/0.018  | +8°    | 9°         |
| UR2 / 0.022/0.018  | +8°    | 9°         |
| UL3 / 0.022/0.018  | 0°     | 11°        |
| UR3 / 0.022/0.018  | 0°     | 11°        |
| UL45 / 0.022/0.018 | -7°    | 0°         |
| UR45 / 0.022/0.018 | -7°    | 0°         |
| LL12 / 0.022/0.018 | -1°    | 0°         |
| LR12 / 0.022/0.018 | -1°    | 0°         |
| LL3 / 0.022/0.018  | -11°   | 5°         |
| LR3 / 0.022/0.018  | -11°   | 5°         |
| LL4 / 0.022/0.018  | -17°   | 0°         |
| LR4 / 0.022/0.018  | -17°   | 0°         |
| LL5 / 0.022/0.018  | -22°   | 0°         |
| LR5 / 0.022/0.018  | -22°   | 0°         |

TOLERANCE ±0.5