




Infinity Plus Brackets



infinity brackets | infinite possibilities



Sandblasted base provides excellent adaptation and adhesion



Features & Benefits

- Infinity Plus low-profile brackets have all the advantages of our Infinity Mini Brackets but also have a mesh base (rather than a cast MIM base) which has also been sandblasted thereby further increasing retention.
- They have a chamfered slot and perfectly round edges for patient comfort.
- Manufactured from the best quality 17-4 Stainless Steel and MIM (Metal Injection Mould) technology.
- Available in mini or standard size.
- Available in .018" / .022" Roth & M.B.T.