


MINI ROTH BRACKETS

in/out			0.68	0.68	0.66	1.06	0.68	0.68	1.06	0.66	0.68	0.68
torque			-7°	-7°	-2°	8°	12°	12°	8°	-2°	-7°	-7°
tip			0°	0°	11°	9°	5°	5°	9°	11°	0°	0°
width			2.8	2.8	3.2	2.8	3.2	3.2	2.8	3.2	2.8	2.8
Tooth	18/22	w/hook	U5RHK	U4RHK	U3RHK	U2R	U1R	U1L	U2L	U3LHK	U4LHK	U5LHK
		w/o hook	U5R	U4R	U3R					U3L	U4L	U5L



in/out			0.68	0.72	0.66	1.06	1.06	1.06	1.06	0.68	0.72	0.68
torque			-22°	-17°	-11°	-1°	-1°	-1°	-1°	-11°	-17°	-22°
tip			0°	0°	5°	0°	0°	0°	0°	5°	0°	0°
width			2.8	2.8	2.8	2.5	2.5	2.5	2.5	2.8	2.8	2.8
Tooth	18/22	w/hook	L5RHK	L4RHK	L3RHK	L2R	L1R	L1L	L2L	L3LHK	L4LHK	L5LHK
		w/o hook	L5R	L4R	L3R					L3L	L4L	L5L

Size	Package	Hooks	
18/22	With laser mark set	3 with Hook	345 with hook

Please note: Minimum quantity order for single brackets is 5 per tooth

