


Ceramic Self-Ligating Brackets



Easy open and shut ensures simple and fast wire changes

Rhodium Coated Clip for Aesthetic effect



Circular Base Treatment: Stable Bonding and Easy Debonding



Ceramic Self-Ligating Brackets

Features & Benefits


- The Nova C Ceramic Self-Ligating (SL) Bracket is made from a durable ceramic material that compares to the strength and performance of metal.
- Superb results are achieved whilst still maintaining enhanced aesthetics.
- The Rhodium Coated Clip and high translucency ensures that the bracket is almost invisible from a distance.
- Low profile bracket ensures patient comfort.
- Available in .022" Roth only.
- Other prescriptions and slot sizes coming soon!

Easy opening

1. Slip Nova C Bracket Opener into the hole on Bracket's Clip.
2. Simply push down Nova C Bracket Opener for opening Clip.



Nova C Opener #1



Nova C Opener #2