

Sapphire Brackets

StarlightTM
SAPPHIRE BRACKETS

Superior transparency
ensures ultimate
aesthetics


Low Profile

Silica Coating on
entire bracket ensures
less breakages

Excellent low friction
capabilities

Rounded edges for
patient comfort

Anatomical base



Superior transparency with rounded edges

Sapphire Brackets

Features & Benefits

- The Starlight Sapphire Bracket is manufactured from monocrystalline alumina.
- This ensures superior clarity and virtually no staining.
- It has a Silica Coating over the whole bracket, not just the slot, which gives excellent sliding characteristics which are comparable to metal brackets. It also gives more durability and fewer breakages and is helpful when de-bonding to avoid splintering.
- Has superior transparency with rounded edges.
- It is low profile for patient comfort, minimal contact to soft tissue and has smooth surface and rounded wings.
- It has a curved base.
- The retentive part of the base has been diminished in size (see picture for circular retentive base) which helps at de-bonding while being sufficient for stable bonding.

Special Base for easy bonding and debonding



Circular Base Treatment:
Stable Bonding and Easy Debonding

Colour guide for easier
bracket positioning

