




© Photo: J Kurtz-Hoffmann

Photo-Mirror-Demister


The Mist-Free Mirror For Dental Photography






FSFDM01 R1495 COD

*Stainless Steel Photo Mirror, Front Surface, double side mirror surface
Occlusal large, Adults*




FSFDM03 R1495 COD

*Stainless Steel Photo Mirror, Front Surface, double side mirror surface
Occlusal standard, Adults*




FSFDM05 R1495 COD

*Stainless Steel Photo Mirror, Front Surface, double side mirror surface
Occlusal small for Children*




FSFDM02 R1495 COD

*Stainless Steel Photo Mirror, Front Surface, double side mirror surface
Lateral, Oral, Adults*



FSFDM04 R1495 COD

*Stainless Steel Photo Mirror, Front Surface, double side mirror surface
Lateral, for Children*



FSFDM06 R1495 COD

*Stainless Steel Photo Mirror, Front Surface, double side mirror surface
Lateral, Universal*

PMDECO R7500 COD

ECONOMY kit consisting of:

- Mirror holder / main unit
- Lithium battery and charger
- Operating instructions
- 1 Photo Mirror of your choice

BULK DEAL

Buy 2 or more
and only pay
R6800 COD



PMDECO+ R8900 COD

ECONOMY PLUS kit consisting of:

- Mirror holder / main unit
- Lithium battery and charger
- Operating instructions
- 1 Photo Mirror of your choice
- 1 Vestibule retractor for CRLA
- 2 Cheek retractors Columbia CRC

BULK DEAL

Buy 2 or more
and only pay
R7900 COD



PMDPRO R21900 COD

PROFESSIONAL kit consisting of:

- Suitcase for photo equipment
- Mirror holder / main unit
- Lithium battery and charger
- Operating instructions
- Photo Mirror Tray PMDK03
- 4 Photo Mirrors of your choice
- 1 Vestibule retractor for CRLA
- 2 Cheek retractors Columbia CRC

BULK DEAL

Buy 2 or more
and only pay
R19750 COD



PMDEVE R16500 COD

EVERY kit consisting of:

- Mirror holder / main unit
- Lithium battery and charger
- Operating instructions
- Photo Mirror Tray PMDK03
- 4 Photo Mirrors of your choice

BULK DEAL

Buy 2 or more
and only pay
R14800 COD



PMDVAR R14900 COD

VARIO kit consisting of:

- Mirror holder / main unit
- Lithium battery and charger
- Operating instructions
- Photo Mirror Tray PMDK03
- 3 Photo Mirrors of your choice

BULK DEAL

Buy 2 or more
and only pay
R13500 COD



PMDVAR+ R16200 COD

VARIO PLUS kit consisting of:

- Mirror holder / main unit
- Lithium battery and charger
- Operating instructions
- Photo Mirror Tray PMDK03
- 3 Photo Mirrors of your choice
- 1 Vestibule retractor for CRLA
- 2 Cheek retractors Columbia CRC

BULK DEAL

Buy 2 or more
and only pay
R14500 COD



MIST-FREE PHOTO MIRRORS

Features and Benefits

Mist free surface of photo mirror due to air stream

- Air volume can be increased or decreased according to patient requirement and/or doctor preference via +/- keypad.
- The air volume level will be indicated on the keypad via LED's.
- The last position of air volume level (after 30 seconds use) will be memorized for the next use.

Angulation of the Photo Demister Handle

- Allows more ergonomic handling.
- Better access for the digital camera.
- Can be easily held by patient (doctor can take photos without an assistant).

Modified mirror clamp fixation system, made of stainless steel

- Secure holding of the mirrors.
- Easy One-Screw-Fixation of the mirror.
- Photo Mirrors can be fixed in the PMD on both sides.

New material of Photo Mirror Demister Cover

- Contains silver ions - antibacterial influence.
- Can be completely disinfected with a disinfectant wipe.
- Will not discolour, even after extended disinfection and use.
- Light weight.
- More robust construction.

The LED light

- Helps with focusing of the digital camera.
- Will not blind patients eyes.
- Light beam completely illuminates oral cavity.

Charging of battery

- Very easy via standard Micro-USB connector.
- Red LED indication of charging progress.
- Battery overload protection.
- Easy battery and charger replacement - available to purchase from Orthoshop.

Stainless steel mirrors

- Better thermal conductivity than glass mirrors – therefore less misting up.
- More robust and longer life than glass mirrors.
- Mirror surface is on both sides.
- Easy to clean and sterilise (also ultrasonic cleaning).

Diverse accessories for more comfort

- Sterilisation cassette for photo mirrors - comfortable sterilisation and storage.
- Retractors made of stainless steel - very long life.

