

LEDEX™

DENTAL CURING LIGHT

WL-090

PERFORMANCE WITH REAL HIGH INTENSITY

-  2 Sec. Cure 2 mm
-  6 Powerful Modes
-  Fast Ortho
-  Turbo 2400 mW/cm²
-  Built-in Radiometer
-  Cartridge Battery



6 powerful extensive modes

OLED display, turbo light guide rod, 6 modes including Low, Ramp, Standard, High, Fast Ortho, Turbo.



Turbo mode - 2 sec. Fast cure

Use high power 10 W LEDs, up to 2400 mW/cm², it can cure the composite resin over 2 mm in 2 seconds.

Never shutdown within continuous use

Smart heat-management system prevents device over-heating during continuous use.

Built-in radiometer

OLED display, easy to check intensity.



Cartridge battery

Use cartridge design, easy to change battery.



Specification

Handpiece		Cradle	Weight	140 g
Light source	High power 10 W LEDs, average lifespan is 100,000 hours		Dimensions	Ø115 x H68 mm
Wavelength	440-480 nm, peak:460 nm	Battery	Type	Lithium-ion, 3.7 V/2500 mAh
Radiant intensity	Up to 2400 mW/cm ²		Charge time	About 2.5 hours
Light guide rod	Ø11>8 mm optical fiber, Steriliz-able in autoclave		Use time	120 Minutes (720 times in 10 Sec. Standard mode)
Dimensions	Ø 38 mm (max.) x L190 mm	Power Supply	Type	Mini USB jack
Weight	180 g (with battery & light guide rod)		Input	100-240 V, 50/60 Hz
			Output	DC 5 V/2 A



Tel: +27 (0)21 976 5403 | Fax: +27 (0)86 538 8694 | Cell: Marius Jooste +27 (0)82 907 9697
Email (Sales): sales@orthoshop.co.za | Skype: orthoshop.sales | WhatsApp: +27 (0)82 907 9697



Made in Taiwan
Patent D124847